



# Digital Microscope

## Instruction Manual



November 2012

718B

# Table of Contents

Important information.....	4
Product description .....	4
Computer requirements .....	5
Windows based PC .....	5
Mac OS based PC .....	5
Package contents .....	6
Product overview.....	7
Product specification .....	8
Software installation .....	9
Windows based PC .....	9
Mac OS based PC .....	9
Getting started.....	10
Assembling the product.....	10
Connecting the device .....	10
Starting the Digital Viewer software.....	11
Windows based PC.....	11
Mac OS based PC .....	11

Focusing.....	12
Using the Digital Viewer software .....	13
System settings menu .....	14
Device setup .....	15
Timed shot setup .....	15
Movie setup .....	15
Save setting.....	15
Language setting.....	16
Advanced settings.....	16
Saved files.....	17
Uninstalling the Digital Viewer software .....	17
Windows based PC .....	18
Mac OS based PC .....	18
Care and maintenance.....	19
Warning .....	19
Legal information.....	21

# Important information

Please read this instruction manual carefully before using this product, and retain these instructions for future reference.

Improvements and changes to this text necessitated by typographical errors, or improvements to the software and/or equipments may be made at any time without notice.

## Product description

This product is a USB powered device which magnifies specimens up to 41x on a 17" computer monitor.

Snapshot of the specimens can be captured using the shutter button located on the top of the device. Video recording is also available.

# Computer requirements

## Windows based PC


Compatible operating systems	Windows XP SP2, SP3 Windows Vista (32 bit or 64 bit) Windows 7 (32 bit or 64 bit) Windows 8 (32 bit or 64 bit)
CPU Speed	P4-1.8GHz or above
RAM	512 MB or above
Hard disk	800 MB or above
USB	USB 2.0

## Mac OS based PC


Compatible operating systems	Mac OS X 10.4.8 - Mac OS X 10.8.1
CPU	Power PC G3/G4/G5 or Intel based
RAM	128 MB or above
Hard disk	800 MB or above
USB	USB 2.0

# Package contents

- |    |                    |    |
|----|--------------------|----|
| 1. | Digital Microscope | x1 |
| 2. | Application CD     | x1 |
| 3. | Instruction manual | x1 |
| 4. | Adaptor (large)    | x1 |
| 5. | Adaptor (small)    | x1 |



# Product overview



- 1. Shutter button
- 2. Focusing ring
- 3. Lens

- 4. LED lights
- 5. Adaptor slot
- 6. USB cable

# Product specification

Connection type	USB 2.0
Effective magnification	14" monitor - 34x 17" monitor - 41x
Effective viewing area	8mm x 6mm
Illumination	Four LED array
Sensor	CMOS
Sensor resolution	680 x 480 pixels (VGA)
Size	75mm x 50mm x 50mm
Weight	133 grams



# Software installation


## Windows based PC

1. Insert the supplied application CD to the CD-ROM of the computer.
2. Double click the "Digital Viewer.exe" icon located within the CD.
3. Follow the Digital Viewer setup wizard to install the software and drivers for the Digital Microscope.




## Mac OS based PC

1. Insert the supplied application CD to the CD-ROM of the computer.
2. Double click the "Digital Viewer.dmg" icon located within the CD.
3. Drag the Digital Viewer icon into the Applications folder.



# Getting started

## Assembling the product



Insert one of the adaptors into the adaptor slot, and fasten it by turning it gently in the clockwise direction.

## Connecting the device

Connect the device to a USB socket of the computer.

## **Starting the Digital Viewer software**

### **Windows based PC**

The Digital Viewer software can be launched by double clicking the Digital Viewer icon from the desktop, or from the start menu (Start > All Programs > Digital Viewer > Digital Viewer).




### **Mac OS based PC**

The Digital Viewer software can be launched by double clicking the Digital Viewer icon from the Applications menu.



## Focusing









Adjust the focus of the image manually by rotating the focusing ring.



# Using the Digital Viewer software

The following buttons are available:

	Clicking this will open the <b>system settings</b> menu (see next section).
	Clicking this will <b>capture</b> the on screen image of the device.
	Clicking this will activate or deactivate the <b>timed shot</b> feature. When activated, images will be captured at a regular interval.
	Clicking this will activate or deactivate the <b>video recording</b> feature.
	Clicking this will show the copyright and version <b>information</b> of the application software.
	Clicking this will <b>shut down</b> the application software.

# System settings menu

The first time the Digital Viewer is opened, the default settings will be loaded, you may change these settings manually in the system settings menu.

The screenshot shows a 'Settings' dialog box with a standard Windows-style title bar (blue with a close button). The dialog is organized into several sections, each with a title and a light gray background:

- Device Setup:** Contains two dropdown menus labeled 'Device:' and 'Resolution:'.
- Timed Shot Setup:** Contains two rows of controls. The first row is for 'Take' (60) shots per minute/hour/day (max. 300). The second row is for 'Duration' (60) minute(s)/hour(s)/day(s) (max. 999). Radio buttons are used to select the time unit.
- Movie Setup:** Contains a 'Resolution:' dropdown menu set to '640X480' and a 'Maximum file size:' text box set to '50' MB (16MB - 4097MB).
- Save Setting:** Contains a text box for the save path ('C:\Users\dneedham\Documents') with a browse button ('...') and a 'Folder:' dropdown menu set to 'Not to create new folder'.
- Language Setting:** Contains a dropdown menu for the language, currently set to 'English'.

On the right side of the dialog, there are four buttons: 'Apply' (highlighted in blue), 'Undo', 'More...', and 'Default'.

## **Device setup**

If the image captured by the Digital Microscope was not displayed by default, you can change this by selecting it from the "Device" drop-down menu.

The resolution of images you capture can be changed from the "Resolution" drop-down menu.

## **Timed shot setup**

The frequency and duration of automatic images capture can be adjusted under this option.

## **Movie setup**

The resolution of videos you record can be changed from the "Resolution" menu. You can also set a maximum file size for each video.

## **Save setting**


The default location for captured images or videos can be changed under this option.

## Language setting

The language of the Digital Viewer software can be changed under this option.

## Advanced settings

By clicking the "More..." button on the right of the system settings menu, you will be able to manually adjust all of the image settings. Note that the settings available may be different, depending on your operating system.





## Saved files

With the Digital Viewer application software opened, you can locate the saved files folder by clicking the "More..." button located on the left of the main software window.



## Full screen

To activate the full screen mode, click the full screen button located on the bottom right hand corner of the Digital Viewer application software window.

To exit full screen mode, either double click on the screen, or press the "Esc" button on the keyboard.

# Uninstalling the Digital Viewer software

## Windows based PC

Select "Uninstall Digital Viewer" from the start menu (Start > All Programs > Digital Viewer > Uninstall Digital Viewer).

## Mac OS based PC

Drag the Digital Viewer application icon from the "Applications" folder to "Trash".

## Care and maintenance

- Keep the device dry and protect it from water or vapour.
- Do not leave your device in a place with an extreme high or low temperature.
- Do not touch the device with a wet hand as it may damage the device, or cause an electric shock to the user.
- Do not use or store the device in dusty, dirty areas as its moving parts may be damaged.
- Do not use harsh chemicals, cleaning solvents or strong detergents to clean the device. Wipe it with a soft cloth slightly dampened in a mild soap-and-water solution.

## Warning

- Do not shine the Digital Microscope's LEDs in the eye, doing so may cause permanent eye damage.

- Do not attempt to open or dismantle the Digital Microscope.

TTS Group Ltd

Park Lane Business Park


Kirkby-in-Ashfield

Nottinghamshire, United Kingdom

NG17 9GU

Freephone: 0800 318686

Freefax: 0800 137525



WARNING! Not suitable for children under 36 months due to small parts. Choking hazard

Product Code: EL00015

Made in China

# Legal information

This document is published without any warranty. While the information provided is believed to be accurate, it may include errors or inaccuracies. In no event shall the manufacturer or its distributors be liable for incidental or consequential damages of any nature, including but not limited to loss of profits or commercial loss, arising out of the use of the information in this document.

Intel is a trademark of Intel Corp. in the U.S. and other countries. Mac, Mac OS and OS X are trademarks of Apple Inc., registered in the U.S. and other countries. PowerPC™ and the PowerPC logo™ are trademarks of International Business Machines Corporation, used under license therefrom. Windows is a registered trademark of Microsoft Corporation in the United States and other countries.